AMERICAN SURGERY

www.americansurgerysa.com

[f /AmericanSurgery](#) [v /AmericanSurgerySA](#) [t /AmericanSurgery](#)

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C-PATCH is a biological matrix intended for the replacement and repair of dura mater. It is obtained from porcine pericardium carefully selected taking into account tissue quality and thickness required. It is subsequently treated with exclusive processes that allow the removal all of the potentially immunogenic elements: cells, lipids, and non-collagenous proteins (NCPs). Our exclusive low temperature production process allows to preserve collagen fibers in their natural conformation. Gamma-ray irradiation guarantees the complete removal of pathogens, without altering the three dimensional structure of the fibers, as well as the preservation of the biomechanical characteristics of the product.

C-Patch works as a bio-inert structure, capable of sustaining tissue repair process without causing adverse reactions in the surrounding tissues. Once applied, C-PATCH acts as a reinforcement and support matrix for fibroblastic infiltration and substrate for the deposit of new collagen. Then it starts to degrade gradually and is replaced by the patient's new tissue.

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If necessary, mold the membrane before hydrating. Hydrate for 10 minutes in sterile saline solution. Apply in the operating room. The membrane can be fixed with sutures.

(Read the instructions for use carefully)

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