
DOMO

ANATOMIC RADIAL HEAD SYSTEM

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RADIAL HEAD REPLACEMENT PROSTHESIS

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DOMO radial head replacement prosthesis is designed to **regain normal function of the patient's elbow and proximal radiocapitellar and radioulnar joints and restore the length of the forearm**, thus allowing for the flexion/extension and elbow and forearm pronosupination.

The prosthesis anatomical design **mimics the shape of the proximal radius and provides optimal elbow joint stability**. All system components are made of titanium alloy pursuant to ASTM F-136-13 Standard.

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Make an incision, 6 to 7 cm long, from the anconeus muscle to the posterior ulna. Proceed to dissect to the joint capsule. The origin of the anconeus muscle can be released subperiosteally and subsequently retracted to allow for adequate exposure of the capsule.

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Resect the head as close to the fracture as possible. The resection limit is 15 mm.

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Initially enter the intramedullary canal with the square tip and then use the rasp.

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Ream the bone surface until at least 60% of the surface is in contact with the reamer.

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Measure the resected head diameter with the diameter gauge. If the measured diameter is between two sizes, select the smaller diameter.

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Insert the trial stem into the medullary canal (TESTER +0). Place the TESTER +3 and TESTER +6 under the plate until the plate reaches the capitellum. It is critical that the coronoid sits flush with the trochlea during the process. The coronoid separated from the trochlea indicates that the collar is too large.

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Insert the stem by turning clockwise until the underside of the implant sits flush with the reamed surface of the radius. Then place the corresponding plate and tap with the plastic tip impactor. Do not damage the polished surface; if there are any scratches or unpolished surfaces, the cartilage will be damaged.

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The set consists of four plates and a wide variety of stems. **The broad range of combinations provides customized solutions featured to address the needs of the most diverse anatomies.**

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Prosthesis indicated for **Comminuted fracture of the radial head requiring resection and Forearm axial instability**

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