**SAI**<sup>®</sup>  
SOUTH AMERICA IMPLANTS

# TÉCNICA QUIRÚRGICA

*SURGICAL TECHNIQUE*

---

Sistema SubTwist BI  
para tenodesis  
*SubTwist BI Tenodesis System*

---

# SUBtwist

## Sistema SubTwist BI para tenodesis

SubTwist BI Tenodesis System

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

*SubTwist BI anchors are designed for fully arthroscopic proximal biceps tenodesis. This system significantly reduces the number of steps and simplifies the biceps tenodesis procedure. Unlike other models of the SubTwist family, the BI anchors are fitted with a PEEK forked tip that is used to direct the biceps tendon to the bottom of the bone hole. The surgeon does not need to pull the tendon out of the body or grasp the tendon with a suture stitch. It is not necessary to know tendon length in advance, as position and tension can be visualized and adjusted prior to final anchor insertion.*

*For drilling, the system features reamers of various diameters with 0.5 mm increments that provide the perfect fit. These reamers are fitted with starter and centering guidewire, thus allowing direct one-step drilling without the need to ream over a guidewire beforehand.*

*The implants are available in two diameters, 5.5mm and 6.5mm. The screws are 10 mm long.*

*The total length of the implant, including the forked tip, is 15 mm.*

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
*Use the tendon fork to shift the biceps tendon to the side. Drill a hole using a centering reamer with a 0.5 to 1.0 mm diameter larger than the selected screw. Drill perpendicular to the bicipital groove at least 20 mm deep. Remove the reamer from the socket manually. Do not use power. Use a shaver to clean soft tissue around the edge of the drill socket to ease screw insertion.*

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
*Use the forked tip of the SubTwist Bi to steer the tendon to the edge of the bone socket.*

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
*Optional: If additional control of the tendon is required, pass the pre-assembled suture into the bifurcated eyelet around the tendon; pass the ends of the suture through the holes in the forked tip. Tense the suture to pull the tendon into the forked eyelet.*

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*Push the tendon to the bottom of the bone socket, while maintaining some light tension on the suture. Verify that the position and tension of the biceps tendon are as desired. Make sure that the anchor body is in full contact with the bone. Hold the thumb pad steady and rotate the driver handle to advance the screw. It is best to stop with the rounded end of the screw slightly protruding. Advance the screw until it is flush with the surface.*

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*Remove and discard the retention sutures and remove the driver. Trim the excess tendon proximally to complete the repair.*

Productos autorizados por la ANMAT PM 2022-13  
*Medical device authorized by ANMAT PM-2022-13*

**DC-120-00**