Técnica Quirúrgica

Surgical Technique

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Distal Biceps Repair

www.floop.com.ar

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
The Floop Biceps Reattachment System consists of a kit designed to facilitate repair of the distal biceps tendon. The kit includes everything you need for a quick and easy repair with optimal results.

The system consists of a Fiber Loop that allows to stitch the tendon quickly and safely, a small button (and its threader) for cortical fixation and a 7 x 10 mm Peek interference screw for optional double fixation. The entire system is sterilized with ethylene oxide.

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
Make a 4 cm transverse incision approximately 3 cm distal to the elbow flexion crease. Identify the lateral antebrachial cutaneous nerve and retract it laterally. Identify the retracted distal end of the biceps tendon and pass it through the incision. Debride the end of the tendon to remove any degenerative diseased tissue. Stitch it with the Fiber Loop suture supplied with the kit.

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
Thread the two ends of the Fiber Loop suture using the threader. Deliver one end of the suture through one side of the biceps button and back through the opposite side again. Thread the other suture end in the same fashion, starting on the opposite side from the first suture end. Make sure the suture ends are not tangled. Pull each suture end simultaneously to ensure that the button slides freely over the sutures.

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
Drill a 3.2 mm bi-cortical tunnel through the radial tuberosity, aiming 30° ulnar to maximize the distance from the posterior interosseous nerve. Use fluoroscopy to confirm drill placement in the radial tuberosity. Then drill a 7/8 mm unicortical tunnel over the 4 mm guide pin and remove the pin with the reamer. Irrigate the wound thoroughly to remove the remains of bone fragments.

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Maintain tension on the sutures and insert the biceps button through both cortices of the radial tuberosity. Closing the clamp will release the button. Pull both ends of the suture to place the button on the radius. Then slowly apply tension to dock the tendon into the bone tunnel. Once the tendon is fully seated, use the Fiber Loop to tighten it and tie a knot. Prepare the 7mm x 10mm interference screw in the driver. For improved safety while blocking, use the threader to deliver the suture through the screw and the screwdriver.

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Insert the screw into the radial side of the bone tunnel. The screw should seat flush with the anterior cortex.

Instrumentos

Instruments

		Qty.
7091-000	Nitinol para interferencial 1,1mm. x 200 mm. Nitinol for interference screw 1.1mm. x 200mm.	1
8741-000	Mecha 4,00 x 150 mm. Sistema biceps Bit 4.00 x 150mm. Biceps System	1
5620-000	Fresa canulada 7,0 mm. Cannulated drill 7.0 mm	1
5621-000	Fresa canulada 8,0 mm. Cannulated drill 8.0 mm	1
5377-000	Atornillador Torx T-25 Canulado Torx T-25 Screwdriver - Cannulated	1
5876-000	Pinza colocadora de Floop Floop placement forceps	1
2043-000	Floop para reinserción del biceps – Kit Floop for biceps reinsertion - Kit	1