

# Advantage<sup>®</sup> Turbo Shaver

A Revolution in Design



Advantage<sup>®</sup> Turbo Shaver



Thorne 5/32" Chuck



Advantage<sup>®</sup> Console



# Advantage® Turbo Plus 10k® Pump

## AN UNBEATABLE COMBINATION

### Advantage Turbo

- 39% increase in torque for better resection performance
- 12,000 rpm for aggressive bone removal
- More ergonomic design to reduce strain
- Bi-Directional blade insertion for operating flexibility
- One-handed collet action to ease blade insertion
- Using Thorne 5/32" Chuck, Turbo can be used to drill shoulder anchor pilot holes



Advantage Turbo



Advantage Console



Thorne 5/32" Chuck

### ORDERING INFORMATION

Description	Cat. No.
Two-Button Advantage Turbo Shaver .....	D9924
No-Button Advantage Turbo Shaver (Footswitch Only) .....	D9920
Thorne 5/32" Chuck.....	D9900
Advantage Console with Cord 115 Volt.....	D3000
Advantage Console with Cord 230 Volt.....	D3000I
Advantage Console with Cord 100 Volt.....	D3000J
10k Console .....	10k
10K Inflow Only Tube Set .....	10k100
10k Inflow/Outflow Tube Set.....	10k150



10k Pump

- Intuitive one-step cassette loading
- Quick setup
- Integrated pole mount works with any standard I.V. pole
- Compatible with video cart

### CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard  
 Largo, FL 33773-4908  
 (727) 392-6464  
 Customer Service: 1-800-237-0169  
 FAX: (727) 399-5256  
 International FAX: +1 (727) 397-4540  
 email: customer\_service@linvatec.com  
 www.linvatec.com

